

Index

Note: Page numbers with prefix ‘*f*’ refers to figures and ‘*t*’ refers to tables.

A

- abdominal organs
 - MRI applications in, 215, 215*f*, 216*f*
- acquisition matrix, 91
- acquisition parameter modification
 - in motion artifact reduction, 126–7, 127*f*, 128*f*
- actively shielded gradient coil, 184
- ADC. *See* analog-to-digital converter (ADC)
- advanced imaging applications, 147–61. *See also specific types, e.g.,* perfusion
 - diffusion, 147–53, 148*f*, 150*f*, 152*f*
 - diffusion-weighted imaging-related, 151–3
 - fMRI, 156–7
 - noble gas imaging, 159–61, 160*f*
 - perfusion, 153–5, 154*f*–6*f*
 - ultra-high field imaging, 158–9
- AET. *See* Applications Entity Title (AET)
- agitated patients
 - MRI in, 221
- aliasing, 105, 107–8, 107*f*–10*f*
- amplitude
 - in RV waveform characterization, 40
- analog-to-digital converter (ADC), 12–13
- Applications Entity Title (AET), 180
- arterial spin labeling (ASL)
 - continuous, 153
 - in perfusion studies, 153
 - pulsed, 153
- artifact(s), 103–25. *See also specific types, e.g.,* motion artifacts
 - “blooming,” 120
 - chemical shift, 108–111, 109*f*, 111*f*, 111*t*
 - coherence, 113, 115–17, 115*f*–18*f*
 - described, 103
 - external, 119–25
 - flow misregistration, 104, 105*f*, 106*f*
 - magnetic susceptibility difference, 117, 118*f*, 120–1, 121*f*
 - motion, 103–5, 105*f*, 106*f*
 - phase cancellation, 110–11, 111*f*, 111*t*
 - radial, 119, 120*f*
 - sequence protocol-related, 105, 107–19
 - truncation, 112, 112*f*–14*f*
- ASL. *See* arterial spin labeling (ASL)

- atom
 - structure of, 2

B

- b* value
 - in diffusion, 148
- base magnetic field
 - safety precautions related to, 197
- benzoxypipronic tetraacetate (BOPTA)
 - in MRI, 191
- black-blood MRA techniques, 135
- BLADE, 134, 134*f*
- “blooming” artifact, 120
- BOLD fMRI studies, 157
- Boltzmann distribution, 7
- BOPTA. *See* benzoxypipronic tetraacetate (BOPTA)
- bowel
 - MRI applications in, 215, 215*f*
- brain
 - MRI applications in, 203, 203*f*, 204*f*
- breast
 - MRI applications in, 209–10
- bright-blood MRA images, 145, 146*f*
- Brownian motion, 147
- bulk magnetic susceptibility, 8–9

C

- Carr–Purcell–Meiboom–Gill (CPMG) technique, 184
- CASL. *See* continuous arterial spin labeling (CASL)
- center frequency
 - in RV waveform characterization, 40
- cerebral ischemia
 - diffusion-weighted imaging in evaluation of, 151–3
- chemical shift
 - in MR, 14–16, 15*f*, 16*f*
 - in MRS, 162–4, 163*f*, 164*f*
- chemical shift artifacts, 108–10, 109*f*
 - of second kind, 110–11, 111*f*, 111*t*
- chemical shift imaging (CSI), 169
- chronic kidney disease (CKD)
 - classification of, 193, 193*t*
- CKD. *See* chronic kidney disease (CKD)
- clinical applications, 200–21
 - abdominal organs, 215, 215*f*, 216*f*

MRI Basic Principles and Applications, Fifth Edition.

Brian M. Dale, Mark A. Brown and Richard C. Semelka.

© 2015 John Wiley & Sons, Ltd. Published 2015 by John Wiley & Sons, Ltd.

clinical applications (*continued*)

brain, 203, 203*f*, 204*f*
 breast, 209–10
 common bile duct, 215, 216*f*
 for ETSE, 220
 examination design considerations in, 202
 general principles of, 200–2
 great vessels, 210–11, 211*f*, 212*f*
 heart, 210–11, 211*f*, 212*f*
 kidneys, 215, 216*f*
 knee, 207*f*, 208*f*
 liver, 211–15, 213*f*, 214*f*
 musculoskeletal, 206, 207*f*–9*f*
 neck, 204–5
 for out-of-phase gradient echo, 219
 pelvis, 215, 217, 217*f*, 218*f*
 in sedated or agitated patients, 221
 shoulder, 209*f*
 for single-echo spin echo sequence, 218–19
 specific sequences and clinical situations, 218–21
 spine, 205–6, 205*f*, 206*f*
 for spoiled gradient echo, 219
 for standard multiecho spin echo, 220
 for STIR imaging, 220
 for *T1*-weighted techniques, 218–20
 fat saturation, 219–20
 for *T2*-weighted techniques, 220–1
 fat saturation, 221
 thorax, 206–9
 coherence artifacts, 113, 115–17, 115*f*–18*f*
 coil(s)
 actively shielded gradient, 184
 phased-array, 58–9, 60*f*, 61*f*
 surface, 59, 60*f*
 common bile duct
 MRI applications in, 215, 216*f*
 composite pulses, 44–5, 45*f*
 computer systems
 in MRI system, 177–80, 178*f*, 187
 continuous arterial spin labeling (CASL), 153
 continuous wave (CW) mode operation
 in RF selective excitation, 39–44, 41*f*–3*f*
 contrast agents, 189–95
 advantages of, 189
 adverse reactions to, 192–3, 193*t*
 BOPTA in, 191
 categorization of, 189
 Feridex, 194–5
 GBCAs, 190–3, 191*f*, 192*t*–3*t*
 intravenous agents, 190–5, 191*f*, 192*t*–3*t*, 194*f*
 negative agents, 195
 oral agents, 195
 safety precautions related to, 199
 T1 relaxation agents, 190–3, 191*f*, 192*t*–3*t*
 T2 relaxation agents, 194–5, 194*f*
 contrast-to-noise ratio, 91
 CPMG technique. *See* Carr–Purcell–Meiboom–Gill (CPMG) technique
 crosstalk, 89–90, 90*f*
 cryogen(s)
 described, 197

safety precautions related to, 197–8

CSI. *See* chemical shift imaging (CSI)

CW. *See* continuous wave (CW)

D

data acquisition system
 in MRI system, 178*f*, 186–8
 data collection methods, 54–5, 54*f*–6*f*
 demodulated signal, 13–14
 DICOM (digital imaging and communications in medicine), 179–80
 DICOM Service Class Uses and Service Class Provider, 180
 diffusion, 147–53, 148*f*, 150*f*, 152*f*
 defined, 147
 described, 147–8, 148*f*
 measurements of, 148–9
 diffusion-weighted imaging
 applications of, 151–3
 digital imaging and communications in medicine (DICOM), 179–80
 directional dependence
 in diffusion measurement, 149–50, 152*f*
 Dixon method
 for fat suppression, 102, 102*f*
 duration
 in RV waveform characterization, 40

E

echo planar imaging (EPI) sequences, 75–7, 76*f*
 MP pulses incorporated into, 82–5, 83*f*–5*f*
 echo spacing, 88
 echo time (*TE*), 87–8, 88*f*
 echo train inversion recovery sequence, 79*f*, 80
 echo train length, 88
 echo train spin echo (ETSE), 69–70, 69*f*
 clinical applications for, 220
 encoding
 phase, 33–5, 33*f*, 34*f*
 energy absorption
 macroscopic, 10, 11*f*
 microscopic, 10, 11*f*
 EPI sequences. *See* echo planar imaging (EPI) sequences
 ETSE. *See* echo train spin echo (ETSE)
 excitation
 RF
 in MR, 10–12, 11*f*
 RF selective, 39–44, 41*f*–3*f*
 excitation angle, 88
 excitation order
 crosstalk and, 89–90, 90*f*
 described, 90
 external artifacts, 119–25
 described, 119
 magnetic field distortions, 120–2, 121*f*, 122*f*
 magnetic susceptibility difference artifact, 120–1, 121*f*
 measurement hardware, 122–4
 noise, 123*f*–4*f*, 124–5

extrinsic parameters, 89–91, 90f

effects of, 92t

variables in, 93

extrinsic variables, 93

F

fat saturation techniques

T1-weighted

clinical applications for, 219–20

T2-weighted

clinical applications for, 221

FDA. *See* Food and Drug Administration (FDA)

Feridex

in MRI, 194–5

ferumoxides, 194–5, 194f

FID. *See* free induction decay (FID)

field of view (FOV), 91

FLAIR. *See* fluid-attenuated inversion recovery (FLAIR)

flip angle, 88

flow compensation

in motion artifact reduction, 132–3, 133f

flow misregistration artifact, 104, 105f, 106f

fluid-attenuated inversion recovery (FLAIR), 81

fMRI. *See* functional MRI (fMRI)

Food and Drug Administration (FDA)

on MRI safety, 196

4D time-resolved MR angiography with keyhole (4D-TRAK), 138

4D-TRAK (4D time-resolved MR angiography with keyhole), 138

Fourier imaging

partial, 53, 53f

Fourier transformation, 103, 169, 171–3, 171f,

172f, 174f, 175f

FOV. *See* field of view (FOV)

free induction decay (FID), 12–13, 13f

frequency encoding, 30–2, 31f, 32f

frequency selective pulses

in RV waveform characterization, 41–4, 42f, 43f

frequency selective saturation, 99–101, 99f, 100f

full-width at half-maximum height (FWHM), 165

functional brain imaging, 156–7

functional MRI (fMRI), 156–7

functional MRI (fMRI) studies

BOLD, 157

FWHM. *See* full-width at half-maximum height (FWHM)

G

gadolinium-based contrast agents (GBCAs)

in MRI, 190–3, 191f, 192t–3t

gating

in motion artifact reduction, 127–32, 130f, 131f

Gaussian pulses

in RV waveform characterization, 42, 43f

GBCAs. *See* gadolinium-based contrast agents (GBCAs)

Generalized Autocalibrating Partially Parallel

Acquisition (GRAPPA), 63–4, 64f

GMR. *See* gradient motion rephasing (GMR)

gradient(s), 182

safety precautions related to, 198

gradient coil

actively shielded, 184

gradient echo sequences, 70–5, 71t, 72f–5f. *See also* specific types

absence of 180° RF pulse in, 70–1, 71t

acronyms for, 71t

described, 70

MP pulses incorporated into, 82–5, 83f–5f

refocused, 71t, 72–5, 72f–5f. *See also* refocused gradient echo sequences

spoiled, 71–2, 71t

steady-state, 71t, 72–5, 72f–5f

types of, 71–5, 71t, 72f–5f

gradient fields, 26–8, 27f, 182

gradient motion rephasing (GMR), 132–3, 133f

gradient pulses, 183

gradient system

in MRI system, 178f, 182–4, 188

GRAPPA (Generalized Autocalibrating Partially Parallel Acquisition), 63–4, 64f

great vessels

MRI applications in, 210–11, 211f, 212f

H

hard pulses, 41, 41f

heart

MRI applications in, 210–11, 211f, 212f

Hermitian symmetry, 53

hydrogen MRS, 162–76. *See also* magnetic resonance spectroscopy (MRS)

hyperbolic secant pulses

in RV waveform characterization, 42, 43f

I

IEC. *See* International Electrotechnical Commission (IEC)

image-based parallel techniques, 62, 63f

image contrast

measurement parameters and, 86–93

image data matrices

raw data in, 46–7, 46f, 47f

image processing computer

in MRI system, 178, 178f

instrumentation, 177–88. *See also* specific components and magnetic resonance imaging (MRI) system

International Electrotechnical Commission (IEC)

on MRI safety, 196

intravenous agents

in MRI, 190–5, 191f, 192t–3t, 194f

intra-voxel incoherent motion, 148–9

intrinsic parameters, 87–9, 87f–9f

effects of, 92t

variables in, 92–3

intrinsic variables, 92–3

inversion recovery (IR) looping modes, 79f, 80

inversion recovery (IR) pulse sequences, 78–80, 78f, 78t

acronyms for, 78t

image reconstruction process in, 80, 81f, 82f

inversion time (TI), 88–9, 89f

IR. *See* inversion recovery (IR)

J*J* coupling

in MRS, 164–5, 166*f*

K*k*-space

raw data and, 48–50, 49*f*, 51*f*

k-space filling techniques, 51–8

reduced, 51–3, 52*f*, 53*f*

reordered, 54–5, 54*f*–6*f*

types of, 56–8, 57*f*–9*f*

k-space formalism, 49–50, 51*f**k*-space trajectory, 54–5, 54*f*–6*f**k*-space-based techniques

GRAPPA, 63–4, 64*f*

parallel, 62, 63*f*

SMASH, 63–4

SNR in, 64

kidney(s)

MRI applications in, 215, 216*f*

knee

MRI applications in, 207*f*, 208*f*

L

Larmor equation, 6, 27

Larmor precession, 4–6, 5*f*

liver

MRI applications in, 211–15, 213*f*, 214*f*

localization techniques

in MRS, 167–9, 168*f*

loop(s)

sequence, 35–8, 36*f*, 37*f*

M

magnet(s)

in MRI system, 178*f*, 180–2, 187

magnet system

in MRI system, 178*f*, 180–2, 187

magnetic field(s), 1–2

magnetic field distortions, 120–2, 121*f*, 122*f*magnetic field gradients, 26–7, 27*f*

magnetic resonance (MR)

chemical shift in, 14–16, 15*f*, 16*f*

concepts of, 10–16

described, 1

RF excitation in, 10–12, 11*f*

RF signal detection in, 12–14, 13*f*, 14*f*

magnetic resonance angiography (MRA), 135–46

approaches to, 135–46

described, 135–6

flow selection in, 135, 136*f*

MIP, 144–6, 145*f*, 146*f*

phase contrast, 141–4, 143*f*–4*f*

problems associated with, 136

techniques, 135–46

time-of-flight, 137–40, 137*f*–42*f*

magnetic resonance imaging (MRI)

advanced imaging applications for, 147–61. *See also specific types and advanced imaging*

applications

in agitated patients, 221

artifacts in images, 103–25. *See also artifact(s)*

clinical applications of, 200–21. *See also clinical applications*

composite pulses in, 44–5, 45*f*

confusing aspects of, 65

contrast agents in, 189–95. *See also contrast agents*

data collection methods, 54–5, 54*f*–6*f*

functional, 156–7

gradient fields in, 26–8, 27*f*

initial step in, 28

k-space filling techniques in, 51–8. *See also*

k-space filling techniques

nonrectilinear sampling in, 56

parallel acquisition methods in, 60–4, 63*f*, 64*f*

phase encoding in, 33–5, 33*f*, 34*f*

phased-array coils in, 58–9, 60*f*, 61*f*

principles of, 26–64

projection reconstruction in, 57, 58*f*

pulse sequences in, 65–85. *See also pulse*

sequence(s)

radial scanning in, 57, 58*f*

ramped sampling in, 56–7, 57*f*

raw data in

image data matrices and, 46–7, 46*f*, 47*f*

k-space and, 48–50, 49*f*, 51*f*

readout or frequency encoding in, 30–2, 31*f*, 32*f*

reduced *k*-space techniques in, 51–3, 52*f*, 53*f*

reordered *k*-space filling techniques in, 54–5, 54*f*–6*f*

RF selective excitation in, 39–44, 41*f*–3*f*

safety precautions in, 196–9. *See also safety precautions*

in sedated patients, 221

sequence looping in, 35–8, 36*f*, 37*f*

slice selection in, 28–9, 28*f*–30*f*

SNR and tradeoffs in, 47–8

two-dimensional spiral sampling in, 58, 59*f*

magnetic resonance imaging (MRI) protocols

principles in development of, 200

magnetic resonance imaging (MRI) system,

177–88. *See also specific components*

actively shielded gradient coil, 184

components list, 187–8

computer systems, 177–80, 178*f*, 187

data acquisition system, 178*f*, 186–8

diagram of, 178*f*

gradient system, 178*f*, 182–4, 188

magnet system, 178*f*, 180–2, 187

RF transmitter system, 178*f*, 184–6, 185*f*, 188

magnetic resonance (MR) signal, 12, 13*f*

magnetic resonance spectroscopy (MRS), 162–76

chemical shift in, 162–4, 163*f*, 164*f*

concepts, 162–7, 163*f*, 164*f*, 166*f*

hydrogen, 162–76

localization techniques, 167–9, 168*f*

multiple voxel techniques, 169, 170*f*

single voxel techniques, 167–9, 168*f*

spectral analysis and postprocessing steps, 169,

171–3, 171*f*, 172*f*, 174*f*, 175*f*

spectral linewidth in, 163*f*, 165–7, 166*f*

spin coupling in, 164–5, 166*f*

SVS in, 167–9, 168*f*

ultra-high field spectroscopy, 173, 175–6, 176f
 water suppression in, 167
 magnetic susceptibility, 8–9
 magnetic susceptibility difference artifact, 117,
 118f, 120–1, 121f
 magnetization materials
 susceptibility and, 8–9
 magnetization-prepared (MP) pulses
 incorporated into gradient echo and EPI
 sequences, 82–5, 83f–5f
 magnetization-prepared (MP) sequences, 77–85,
 78f, 78t, 79f, 81f–4f
 example of, 78, 78f, 78t
 magnetization transfer (MT) pulse, 96
 in time-of-flight MRA, 138
 magnetization transfer (MT) suppression, 96–9,
 97f, 98f
 MAST. *See* motion artifact suppression techniques
 (MAST)
 maximum intensity projection (MIP)
 MRA-related, 144–6, 145f, 146f
 maximum signal content
 site of, 48
 measurement controller
 in MRI system, 178, 178f
 measurement hardware
 as artifact, 122–4
 measurement parameters, 86–93
 extrinsic parameters, 89–91, 90f
 effects of, 92t
 intrinsic parameters, 87–9, 87f–9f
 effects of, 92t
 modification of, 86
 parameter tradeoffs, 91–3, 92t
 measurement protocol
 criteria in determining, 91–2
 measurement techniques
 differentiating, 35–7
 MIP. *See* maximum intensity projection (MIP)
 modulation
 in RV waveform characterization, 40
 motion
 Brownian, 147
 intravoxel incoherent, 148–9
 motion artifact(s), 103–5, 105f, 106f
 appearance of, 104
 causes of, 103–4
 flow misregistration artifact, 104, 105f, 106f
 respiratory, 104–5, 106f
 types of, 104, 105f, 106f
 motion artifact reduction techniques, 126–34
 acquisition parameter modification, 126–7, 127f,
 128f
 flow compensation, 132–3, 133f
 radial-based motion compensation, 134, 134f
 triggering gating, 127–32, 130f, 131f
 motion artifact suppression techniques (MAST)
 in motion artifact reduction, 132–3, 133f
 MP sequences. *See* magnetization-prepared (MP)
 sequences
 MR. *See* magnetic resonance (MR)
 MRA. *See* magnetic resonance angiography (MRA)

MRI. *See* magnetic resonance imaging (MRI)
 MRS. *See* magnetic resonance spectroscopy (MRS)
 MT. *See* magnetization transfer (MT)
 multiecho spin echo sequences, 67, 68f
 multiple voxel techniques
 in MRS, 169, 170f
 musculoskeletal system
 MRI applications in, 206, 207f–9f
N
 nature of synthesizer mixing
 in RV waveform characterization, 40
 neck
 safety precautions related to, 204–5
 negative agents
 in MRI, 195
 net magnetization
 described, 6–7, 7f, 8f
 production of, 1–9
 NMR. *See* nuclear magnetic resonance (NMR)
 noble gas imaging, 159–61, 160f
 noise
 as artifact, 123f–4f, 124–5
 nonrectilinear sampling, 56
 nonsaturation methods, 101–2, 101f, 102f
 nonselective pulses
 in RV waveform characterization, 41, 41f
 nuclear magnetic moments, 3f, 4
 nuclear magnetic resonance (NMR), 1
 nuclear spin, 2–4, 3f, 3t, 5f
 nucleus(i)
 spin of, 2–4, 3f, 3t
 null time for tissue, 80
 Nyquist frequency, 14, 32

O
 oral agents
 in MRI, 195
 out-of-phase gradient echo
 clinical applications for, 219
 oversampling, 91
 effects of, 107, 107f

P
 PACS (picture archiving and communications
 system) unit, 178–9
 parallel acquisition methods, 60–4, 63f, 64f
 requirements for, 62, 64f
 types of, 62, 63f
 partial Fournier imaging, 53, 53f
 pelvis
 MRI applications in, 215, 217, 217f, 218f
 Perflubron (perfluorocetyl bromide [PFOB]), 195
 perfluorocetyl bromide (PFOB), 195
 perfusion, 153–5, 154f–6f
 perfusion blood flow
 in diffusion measurement, 148–9
 perfusion studies
 applications of, 154–5, 154f–6f
 approaches to, 153
 PFOB. *See* perfluorocetyl bromide (PFOB)
 PGSE. *See* pulsed gradient spin echo (PGSE)
 method

phase
 in RV waveform characterization, 40
 phase cancellation artifact, 110–11, 111f, 111t
 phase contrast MRA, 141–4, 143f–4f
 phase encoding
 in MRI, 33–5, 33f, 34f
 phased-array coils, 58–9, 60f, 61f
 picture archiving and communications system (PACS) unit, 178–9
 point resolved spectroscopy (PRESS), 167–9, 168f
 postprocessing
 in MRS, 169, 171–3, 171f, 172f, 174f, 175f
 pregnancy
 MRI during, 196
 PRESS. *See* point resolved spectroscopy (PRESS)
 projection reconstruction, 57, 58f
 PROPELLER, 134, 134f
 PSAL. *See* pulsed arterial spin labeling (PSAL)
 pseudodiffusion, 149
 pulse(s)
 composite, 44–5, 45f
 frequency selective, 41–4, 42f, 43f, 99–101, 99f, 100f
 frequency selective soft, 41–4, 42f, 43f
 Gaussian, 42, 43f
 gradient, 183
 hard, 41, 41f
 hyperbolic secant, 42, 43f
 MP
 incorporated into gradient echo and EPI sequences, 82–5, 83f–5f
 MT, 96–9, 97f, 98f
 in time-of-flight MRA, 138
 nonselective, 41, 41f
 in RV waveform characterization, 41, 41f
 RF excitation, 44
 in RV waveform characterization, 39–44, 41f–3f.
 See also specific types
 sinc, 41–2, 42f, 43f
 soft, 41–4, 42f, 43f
 spatial presaturation, 94–6, 95f
 pulse sequence(s), 65–85. *See also specific types*
 comparison of, 65–7, 66f
 defined, 65
 EPI sequences, 75–7, 76f
 gradient echo sequences, 70–5, 71t, 72f–5f
 implementation of, 65
 IR pulse sequences, 78–80, 78f, 78t
 MP sequences, 77–85, 78f, 78t, 79f, 81f–4f
 spin echo sequences, 67–70, 67t, 68f, 69f
 timing diagrams in comparison of, 66–7, 66f
 pulse shape
 in RV waveform characterization, 41, 41f
 pulsed arterial spin labeling (PSAL), 153
 pulsed gradient spin echo (PGSE) method, 148, 148f

R

radial artifact, 119, 120f
 radial-based motion compensation
 in motion artifact reduction, 134, 134f
 radial scanning, 57, 58f

radiofrequency (RF) excitation
 in MR, 10–12, 11f
 selective, 39–44, 41f–3f
 radiofrequency (RF) excitation pulses
 criteria for, 44
 radiofrequency (RF) power deposition
 safety precautions related to, 198–9
 radiofrequency (RF) signal detection
 in MR, 12–14, 13f, 14f
 radiofrequency (RF) transmitter system
 in MRI system, 178f, 184–6, 185f, 188
 radiofrequency (RF) waveform
 parameters in characterizing, 39–44, 41f–3f
 ramped sampling, 56–7, 57f
 raw data, 46–7, 46f, 47f
 k-space and, 48–50, 49f, 51f
 readout, 30–2, 31f, 32f
 reduced *k*-space techniques, 51–3, 52f, 53f
 refocused gradient echo sequences, 71t, 72–5, 72f–5f
 described, 72
 signals measured from, 73–5, 73f, 74f
 relaxation, 17–25
 defined, 17
 spin–spin, 21–2, 21f
 *T*₁
 saturation and, 17–21, 18f–20f
 *T*₂, 21–5, 21f, 22f, 25f
 *T*₂*, 22f, 23–4, 25f
 relaxation agents
 in MRI. *See specific agents and* contrast agents
 reordered *k*-space filling techniques, 54–5, 54f–6f
 repetition time (*TR*), 87, 87f
 respiratory motion artifact, 104–5, 106f
 RF. *See* radiofrequency (RF)
 rise time, 183

S

safety precautions, 196–9
 base magnetic field, 197
 contrast agents, 199
 cryogenics, 197–8
 FDA on, 196
 gradients, 198
 IEC on, 196
 pregnancy-related, 196
 RF power deposition, 198–9
 saturation
 *T*₁ relaxation and, 17–21, 18f–20f
 scanner controller
 in MRI system, 177–8, 178f
 “schimming,” 122
 sedated patients
 MRI in, 221
 “seed” region, 146
 Send/Receive, 180
 sequence looping, 35–8, 36f, 37f
 sequence protocol-related artifacts, 105, 107–19
 aliasing, 105, 107–8, 107f–10f
 chemical shift artifacts, 108–10, 109f
 coherence artifacts, 113, 115–17, 115f–18f

- magnetic susceptibility difference artifact, 117, 118f
 - phase cancellation artifact, 110–11, 111f, 111t
 - radial artifact, 119, 120f
 - truncation artifacts, 112, 112f–14f
 - shoulder
 - MRI applications in, 209f
 - signal suppression techniques, 94–102
 - Dixon method, 102, 102f
 - frequency selective saturation, 99–101, 99f, 100f
 - MT suppression, 96–9, 97f, 98f
 - nonsaturation methods, 101–2, 101f, 102f
 - spatial presaturation, 94–6, 95f
 - water excitation, 101, 101f
 - signal-to-noise ratio (SNR), 14
 - increase in, 91
 - in *k*-space–based techniques, 64
 - in parallel acquisition methods, 61–2
 - tradeoffs, 47–8
 - Simultaneous Acquisition of Spatial Harmonics (SMASH), 63–4
 - sinc pulses
 - in RV waveform characterization, 41–2, 42f, 43f
 - single-echo spin echo sequence, 67, 68f
 - clinical applications for, 218–19
 - single voxel spectroscopy (SVS)
 - in MRS, 167–9, 168f
 - single voxel techniques
 - in MRS, 167–9, 168f
 - slice gap, 89–90, 90f
 - slice selection
 - in MRI, 28–9, 28f–30f
 - slice thickness (*TH*), 89
 - SMASH (Simultaneous Acquisition of Spatial Harmonics), 63–4
 - SNR. *See* signal-to-noise ratio (SNR)
 - soft pulses
 - in RV waveform characterization, 41–4, 42f, 43f
 - spatial presaturation, 94–6, 95f
 - spatial presaturation pulse, 94–6, 95f
 - spectral analysis
 - in MRS, 169, 171–3, 171f, 172f, 174f, 175f
 - spectral linewidth
 - in MRS, 163f, 165–7, 166f
 - spike(s)
 - defined, 123f, 125
 - spin coupling
 - in MRS, 164–5, 166f
 - spin echo(es), 21–5, 21f, 22f, 25f
 - spin echo sequences, 67–70, 67t, 68f, 69f. *See also specific types*
 - acronyms for, 67t
 - ETSE, 69–70, 69f
 - multiecho, 67, 68f
 - single-echo, 67, 68f
 - types of, 67–8, 67t
 - spin echo–based sequences, 82–3
 - spine
 - MRI applications in, 205–6, 205f, 206f
 - spin–spin relaxation, 21–2, 21f
 - SPIO. *See* superparamagnetic iron oxide (SPIO)
 - spoiled gradient echo
 - clinical applications for, 219
 - spoiled gradient echo sequences, 71–2, 71t
 - standard multiecho spin echo
 - clinical applications for, 220
 - steady-state gradient echo sequences, 71t, 72–5, 72f–5f. *See also* refocused gradient echo sequences
 - STEAM. *See* stimulated echo acquisition method (STEAM)
 - Stejskal–Tanner method, 148, 148f
 - stimulated echo acquisition method (STEAM), 167–8, 168f
 - STIR imaging
 - clinical applications for, 220
 - stroke
 - diffusion-weighted imaging in evaluation of, 151–3
 - subloop(s)
 - in sequence looping, 35–6
 - superparamagnetic iron oxide (SPIO) particulate molecules, 194–5, 194f
 - surface coils, 59, 60f
 - susceptibility
 - bulk magnetic, 8–9
 - magnetic, 8–9
 - magnetization materials and, 8–9
 - SVS. *See* single voxel spectroscopy (SVS)
- T**
- TCP/IP (Transmission Control Protocol/Internet Protocol)
 - in MRI system, 179
 - TE*. *See* echo time (*TE*)
 - TH*. *See* slice thickness (*TH*)
 - thorax
 - MRI applications in, 206–9
 - TI*. *See* inversion time (*TI*)
 - time-of-flight MRA, 137–40, 137f–42f
 - time-resolved angiography with interleaved stochastic trajectories (TWIST), 138
 - time-resolved imaging of contrast kinetics (TRICKS), 138
 - timing diagram
 - in pulse sequence comparison, 66–7, 66f
 - T1* relaxation
 - defined, 17
 - saturation and, 17–21, 18f–20f
 - T1* relaxation agents
 - in MRI, 190–3, 191f, 192t–3t
 - T1*-weighted imaging studies
 - in perfusion studies, 153
 - T1*-weighted techniques
 - clinical applications for, 218–20. *See also specific techniques*
 - fat saturation
 - clinical applications for, 219–20
 - TR*. *See* repetition time (*TR*)
 - Transmission Control Protocol/Internet Protocol (TCP/IP)
 - in MRI system, 179
 - TRICKS (time-resolved imaging of contrast kinetics), 138

triggering gating
 in motion artifact reduction, 127–32, 130*f*,
 131*f*
truncation artifacts, 112, 112*f*–14*f*
T2 effects
 in diffusion measurement, 149, 150*f*
T2 relaxation, 21–5, 21*f*, 22*f*, 25*f*
*T2** relaxation, 22*f*, 23–4, 25*f*
 defined, 23
T2 relaxation agents
 in MRI, 194–5, 194*f*
T2-weighted imaging studies
 in perfusion studies, 153
T2-weighted techniques
 clinical applications for, 220–1
 fat saturation
 clinical applications for, 221
turbo factor, 88

TWIST (time-resolved angiography with
 interleaved stochastic trajectories), 138
two-dimensional spiral sampling, 58, 59*f*

U

ultra-high field imaging, 158–9
ultra-high field spectroscopy
 in MRS, 173, 175–6, 176*f*

W

water excitation
 for signal suppression, 101, 101*f*
water suppression
 in MRS, 167

Z

Zeeman interaction, 6–7, 7*f*
zero filling, 171, 171*f*