

Glossary

See also: Liney G (2011) *MRI from A to Z*, 2nd edition. London: Springer-Verlag.

Acronym/symbol	Term	Explanation
γ	Gyromagnetic ratio	A constant which relates the precessional frequency of a nucleus to the external magnetic field B_0 . For ^1H $\gamma = 2.68 \times 10^8$ radians s^{-1} .
μ	Magnetic moment	Microscopic magnetic field originating from nuclear spin.
α	Flip angle	The angle through which magnetization M is tipped, or flipped, by an RF B_1 pulse.
ω_0, f_0	Larmor frequency	The precession frequency of the nuclear spins, or the resonance frequency for nuclear magnetic transitions.
ΔB	Inhomogeneity	Non-uniformity (of a magnetic field), expressed in parts-per-million (ppm).
τ_c	Correlation time	The average time between molecular collisions, part of the relaxation mechanism of protons.
χ_m	Susceptibility	Property determining the magnetization of a material in an external magnetic field.
180°	180° pulse	An RF (radiofrequency) pulse which flips the magnetization through 180° . May be used either for inversion or refocusing.
2D FT	Two-Dimensional Fourier Transform	The process whereby the data from frequency and phase-encoded MR signals is converted to a two-dimensional image.
3D FT	Three-Dimensional Fourier Transform	A volume-based Cartesian image acquisition which utilizes two phase-encode directions.
4D-TRAK	4D-Time-Resolved Angiography using Keyhole	Method for time-resolved CE-MRA (Philips)
90°	90° pulse	An RF (radiofrequency) pulse which flips the magnetization through 90° . Usually used for excitation, but can also act as a refocusing pulse.
ARC	Auto-calibrating Reconstruction for Cartesian sampling	A k-space-based parallel imaging method that uses auto-calibration, related to GRAPPA (GE Healthcare).
ASL	Arterial Spin Labelling	A method to assess tissue perfusion using endogenous water as a contrast agent.

(cont.)

Acronym/symbol	Term	Explanation
ASSET	Array Spatial Sensitivity Encoding Technique	An image-based parallel imaging method related to SENSE (GE Healthcare).
B	Magnetic field	Magnetic flux density or induction, measured in tesla (T).
B₀		The main (static) magnetic field, e.g. 1.5 T.
B₁		The alternating radiofrequency magnetic field used to rotate the net magnetization.
BASG	Balanced SARGE	A fully rewound GE imaging sequence (Hitachi).
bFFE	Balanced FFE	A fully rewound GE imaging sequence (Philips).
BLADE	BLADE	A hybrid Cartesian/radial acquisition method that can partly correct for patient motion (Siemens; 'BLADE' is not an acronym, it is a tradename). See also PROPELLER.
BLAST	Broad-use Linear Acquisition Speed-up Technique	A method for the acceleration of dynamic images (Philips).
BOLD	Blood Oxygen Level Dependent	Effect of deoxygenated blood on MR signals.
BRAVO	BRAin VOlume imaging	IR-prepared fast 3D gradient echo for isotropic brain imaging (GE Healthcare).
BSI	Blood Sensitive Imaging	Method for non-contrast MR angiography based on 3D TSE (Hitachi).
bSSFP	Balanced Steady-State-Free-Precession	An alternative generic name for fully rewound GE sequences.
b-TRANCE	Balanced Triggered Angiography Non Contrast Enhanced	Method for non-contrast MR angiography based on fully rewound GE (Philips).
BW, RBW	Bandwidth, receive bandwidth	The range of frequencies contained in a pulse or signal. The range of frequencies that may be detected by the MR receiver.
CARE Bolus	Combined Applications to Reduce Exposure Bolus	Method for fluoro-triggered contrast-enhanced MR angiography (Siemens).
CE-MRA	Contrast Enhanced Magnetic Resonance Angiography	Methods for acquiring MR angiograms using an injection of gadolinium-based contrast agent.
CENTRA	Contrast ENhanced Timing Robust Angiography	Centric k-space method for time-resolved CE-MRA (Philips).
CHESS	CHEmical Shift Selective	Chemical excitation or saturation of either water or fat ¹ H nuclei based upon differences in the Larmor frequency of each chemical species.
CIA	Contrast Improved Angiography	Method for non-contrast MR angiography based on 3D TSE (Toshiba).
CISS	Constructive Interference in the Steady State	A fully rewound GE imaging sequence with two acquisitions and phase-cycling to reduce banding artefacts (Siemens).

(cont.)

Acronym/symbol	Term	Explanation
CLEAR	Constant LLevel AppeaRance	Post-processing filter to reduce signal inhomogeneities caused by coil sensitivity (Philips).
CP	Carr–Purcell sequence	Sequence containing a series of 90° – 180° – 180° –... pulses, resulting in an echo train of length ETL.
CPMG	Carr–Purcell–Meiboom–Gill sequence	Similar to CP but with a 90° phase difference between the excitation and refocusing pulses. Even-numbered echoes are corrected for inaccuracies in the refocusing pulses.
CS	Compressed sensing	Rapid image acquisition technique applicable to sparse data sets.
CUBE		3D TSE sequence with variable refocusing flip angles (GE Healthcare; Cube is not an acronym, it is a tradename).
DCE	Dynamic Contrast Enhanced	A dynamic T_1 w method used to monitor the signal changes as a contrast agent passes through tissue.
DESS	Double Echo Steady State	A form of GE imaging sequence that combines two different types of gradient-echo signals to form one image (Siemens).
DIET	Delayed Interval Echo Train	A TSE imaging sequence with an echo spacing designed to maintain J-coupling, thus with reduced fat signal compared to other TSE sequences (Toshiba).
DIR	Double Inversion Recovery	Method to suppress the signal from flowing blood, also known as black or dark-blood preparation. Alternatively, a dual inversion prepulse used with 3D TSE to null two tissues, typically CSF and grey matter in brain imaging.
Dixon		TSE or GE techniques for producing separate water and fat images, named after the inventor.
DRIVE	DRIVen Equilibrium	90° RF pulse added to end of FSE echo train to restore longitudinal magnetization (Philips).
DRKS	Differential Rate K-space Sampling	Centric k-space method for time-resolved CE-MRA (Toshiba).
DSC	Dynamic Susceptibility Contrast	A dynamic T_2^* w method used to monitor the signal changes as a contrast agent passes through tissue.
DSV	Diameter Spherical Volume	Spherical region with this diameter over which B_0 uniformity is measured.
DTI	Diffusion Tensor Imaging	Method to determine the anisotropy of tissue diffusion.
EPI	Echo Planar Imaging	A pulse sequence that allows rapid MR acquisition, where k-space is collected following a single excitation pulse.
ESP	Echo SPacing	The time interval between refocusing pulses in a TSE sequence.

(cont.)

Acronym/symbol	Term	Explanation
ETL	Echo Train Length	The number of echoes acquired following a single initial excitation.
FASE	Fast Advanced Spin Echo	A multiple spin-echo sequence that uses a segmented k-space (Toshiba). SuperFASE is a variant with very short echo spacing. See also TSE.
FBI	Fresh Blood Imaging	Method for non-contrast MR angiography based on 3D TSE (Toshiba).
FE	Frequency Encoding	A process using a frequency encode gradient to determine the position of MR signal during its acquisition.
FFE, TFE	Fast Field Echo, Turbo Field Echo	A GE imaging sequence (Philips). See also T1-FFE, bFFE, T2-FFE. TFE is a fast version of FFE.
FID	Free Induction Decay	MR signal following a single RF excitation pulse.
FIESTA	Fast Imaging with Enhanced Steady sTate Acquisition	A fully rewound GE imaging sequence (GE Healthcare).
FIESTA-c	Fast Imaging with Enhanced Steady sTate Acquisition Cycled phases	A fully rewound GE imaging sequence with two acquisitions and phase-cycling to reduce banding artefacts (GE Healthcare). See also CISS.
FIR	Fast Inversion Recovery	A TSE imaging sequence with an inversion prepulse, used for FLAIR or STIR (Hitachi).
FISP, TrueFISP	Fast Imaging with Steady Precession	A rewound GE imaging sequence (Siemens). The fully rewound version is called TrueFISP.
FLAIR	Fluid Attenuated Inversion Recovery	An SE-based imaging sequence used to null the signal from cerebrospinal fluid (CSF) based upon its longitudinal relaxation time T_1 .
FLASH	Fast Low Angle SHot	A spoiled GE imaging sequence using RF spoiling (Siemens). Turbo-FLASH is a fast version of FLASH.
FLUTE	FLUoro TriggEred	Method for fluoro-triggered contrast-enhanced MR angiography (Hitachi).
fMRI	Functional Magnetic Resonance Imaging	Method to investigate signal changes based upon the BOLD contrast mechanism.
f_N	Nyquist frequency	The minimum frequency required to accurately sample (digitize) a signal.
FOV	Field Of View	The size of the region to be imaged.
FSBB	Flow Sensitive Black Blood	Method for susceptibility-weighted imaging (Toshiba).
FSE	Fast Spin Echo	A multiple spin-echo sequence that uses a segmented k-space. See also TSE.
FWHM	Full Width at Half Maximum	A measure of slice width or thickness.
GE	Gradient Echo	A pulse sequence which creates a signal at time TE by refocusing using gradients. Alternatively, the signal itself at time TE.

(cont.)

Acronym/symbol	Term	Explanation
GMN	Gradient Moment Nulling	A method to eliminate the phase shifts due to moving spins. Also known as flow or velocity compensation.
GRAPPA	GeneRALized Auto-calibrating Partial Parallel Acquisition	A k-space-based parallel imaging method that uses auto-calibration (Siemens).
GRASE	Gradient and Spin Echo	A pulse sequence that acquires both spin and gradient echoes within a single TSE-type readout.
GRE	Gradient Recalled Echo	A rewind GE imaging sequence (GE Healthcare).
G_x, G_y, G_z	Gradients	Magnetic field gradients, or linear variations in B_z along orthogonal directions x, y and z , measured in mT m^{-1} .
HASTE	Half-Fourier Acquired Single shot Turbo spin Echo	A single shot TSE imaging sequence that also uses half-Fourier data acquisition and reconstruction. See also SS-TSE.
IDEAL	Iterative Decomposition of water and fat using Echo Asymmetry and Least-squares estimation	A technique for producing separate water and fat images (GE Healthcare). See also Dixon.
IFIR	Inhance inFlow Inversion Recovery	Method for non-contrast MR angiography based on fully rewind GE (GE Healthcare).
iPAT	Integrated Parallel Imaging Techniques	An alternative name for parallel imaging techniques (Siemens).
IR	Inversion Recovery	A pulse sequence that initially inverts M_0 and measures its longitudinal relaxation after inversion time T_I .
ISCE	Inclined Slab for Contrast Enhancement	Linear flip angle variation used in 3D-MRA to reduce saturation effects in slice-encoding direction (Toshiba). See also TONE.
JET		A hybrid Cartesian/radial acquisition method that can partly correct for patient motion (Toshiba; JET is not an acronym, it is a tradename). See also PROPELLER.
LAVA	Liver Acquisition with Volume Acceleration	3D T_1 w GE with fat suppression used for dynamic liver imaging (GE Healthcare). See also THRIVE, VIBE.
M	Magnetization	The magnetic field vector produced in a material when placed in an external magnetic field, with units amperes per meter (A m^{-1}).
M	Gradient moment	The time integral of a gradient waveform.
M_0	Equilibrium Magnetization	Equilibrium magnetization formed from the vector sum nuclear magnetic moments when placed in an external magnetic field B_0 .
MAVRIC	Multi-Acquisition with Variable Resonances Image Combination	Method based on 3D TSE to reduce artefacts when imaging in the presence of metal implants (GE Healthcare). Also a generic term.

(cont.)

Acronym/symbol	Term	Explanation
MEDIC	Multi-Echo Data-Image Combination	Method that combines multiple gradient echoes to improve image quality and contrast (Siemens).
MERGE	Multi-Echo Recombined Gradient Echo	Method that combines multiple gradient echoes to improve image quality and contrast (GE Healthcare).
mFFE	Multi-echo FFE	Method that combines multiple gradient echoes to improve image quality and contrast (Philips).
MIP	Maximum Intensity Projection	An image processing operation used in MR angiography to produce 2D projection images from 3D raw image data.
mIP	minimum Intensity Projection	An image processing operation used in susceptibility-weighted angiography to produce 2D projection images from 3D raw image data.
MOTSA	Multiple Overlapping Thin Slab Acquisition	A technique for acquiring multiple 3D slabs for 3D non-contrast MRA, avoiding saturation effects.
MPR	Multi-Planar Reformat	An image processing operation used to produce new 2D slices at arbitrary orientations from 3D raw image data.
MP-RAGE	Magnetization-Prepared Rapid Acquisition by Gradient Echoes	IR-prepared fast 3D gradient echo for isotropic brain imaging (Siemens). Also a generic term.
MRA	MR Angiography	Imaging of blood vessels, usually arteries.
mSENSE	Modified SENSitivity Encoding	An image-based parallel imaging method related to SENSE using auto-calibration (Siemens).
MT	Magnetization Transfer	Excitation of the bound fraction of nuclei using an off-resonance B_1 pulse.
MultiVane	MultiVane	A hybrid Cartesian/radial acquisition method that can partly correct for patient motion (Philips; MultiVane is not an acronym, it is a tradename). See also PROPELLER.
M_{xy}	Transverse magnetization	Component of M in the xy plane. M_{xy} produces the signal.
M_z	Longitudinal magnetization	Component of M in the z axis, static field direction, denoted M_z .
NATIVE	Non-contrast MR of ArTerles and VEins	Methods for non-contrast MR angiography (Siemens). NATIVE-SPACE is based on 3D TSE; NATIVE-TrueFISP is based on fully rewound GE.
NATURAL	NATural Uniformity Realization ALgorithm	Post-processing filter to reduce signal inhomogeneities caused by coil sensitivity (Hitachi).
NC-MRA	Non-Contrast Magnetic Resonance Angiography	Methods for acquiring MR angiograms without using exogenous contrast agent.
NEX	Number of EXcitations	Number of signal averages (GE Healthcare). See also NSA.
NMR	Nuclear Magnetic Resonance	The basic phenomenon underlying MRI.

(cont.)

Acronym/symbol	Term	Explanation
NSA	Number of Signal Acquisitions	Number of signal acquisitions averaged to improve signal-to-noise ratio (SNR). See also NEX.
O-MAR	Ortho-Metal Artefact Reduction	Method based on multi-slice TSE to reduce artefacts when imaging in the presence of metal implants (Philips). See also SEMAC.
PACE	Prospective Acquisition CorrEction	Technique for updating scan geometry in real time based on detecting motion by navigators (Siemens).
PBSG	Phase Balanced SARGE	A fully rewound GE imaging sequence with two acquisitions and phase-cycling to reduce banding artefacts (Hitachi).
PC	Phase contrast	A non-contrast MRA technique that relies on the phase shift caused by spins moving through a magnetic field gradient.
PD	Proton density	The MR signal intensity in the absence of relaxation, related to number of hydrogen nuclei per unit volume, equivalent to M_0 .
PE	Phase Encoding	A process using a phase-encode gradient to encode the MR signal in terms of spatial frequencies.
PEAKS	PEak Arterial K Space	Centric k-space method for time-resolved CE-MRA (Hitachi).
PEAR	Phase-Encoded Artefact Reduction	Reordered phase-encoding to reduce respiratory motion artefacts (Philips). See also ROPE.
PERRM	Phase-Encode Reordering to Reduce Motion	Reordered phase encoding to reduce respiratory motion artefacts (Hitachi). See also ROPE.
PNS	Peripheral nerve stimulation	Peripheral nerve stimulation is a bio-effect caused by rapidly time-varying magnetic fields.
ppm	Parts Per Million	A measure of the chemical shift between metabolites. Alternatively a measure of field inhomogeneity.
Promo	Prospective Motion correction	Technique for updating scan geometry in real time based on detecting motion by navigators (GE Healthcare).
PROPELLER	Periodically Rotated Overlapping Parallel Lines with Enhanced Reconstruction	A hybrid Cartesian/radial acquisition method that can partly correct for patient motion (GE Healthcare). Also a generic term.
PSIF		A time-reversed GE imaging sequence with improved T_2w contrast (Siemens; PSIF is not an acronym, it is 'FISP' backwards).
PURE	Phased array UnifoRmity Enhancement	Post-processing filter to reduce signal inhomogeneities caused by coil sensitivity (GE Healthcare).
r_1, r_2	Relaxivity	The relaxivity of a contrast agent determines its shortening effect on relaxation times.

(cont.)

Acronym/symbol	Term	Explanation
R_1 , R_2 or R_2^*	Relaxation rate	Reciprocal of relaxation time.
RADAR	RADial Acquisition Regime	A hybrid Cartesian/radial acquisition method that can partly correct for patient motion (Hitachi). See also PROPELLER.
RAPID, RAPID-3D	Rapid Acquisition through Parallel Imaging Design	An image-based parallel imaging method related to SENSE (Hitachi).
rb	Rheobase	The minimum threshold for a physiological response to a stimulus.
RESTORE		90° RF pulse added to end of FSE echo train to restore longitudinal magnetization (Siemens; RESTORE is not an acronym, it is a tradename).
ROPE	Respiratory-Ordered Phase Encoding	Reordered phase encoding to reduce respiratory motion artefacts.
RSSG	RF Spoiled SARGE	A spoiled GE imaging sequence using RF spoiling (Hitachi).
SAR	Specific Absorption Rate	The RF power per unit mass deposited in tissue in $W\ kg^{-1}$.
SARGE	Steady state Acquisition Rewound Gradient Echo	A rewound GE imaging sequence (Hitachi).
SE	Spin echo	A pulse sequence which creates a signal at time TE by refocusing from two RF pulses, usually a 90° followed by a 180°. Alternatively, the signal itself at time TE.
SEMAC	Slice Encoding for Metal Artefact Correction	Method based on multi-slice TSE to reduce artefacts when imaging in the presence of metal implants.
SENSE	SENSitivity Encoding	Technique for parallel imaging based upon image-space unwrapping (Philips). Also a generic term.
SMASH	SiMultaneous Acquisition of Spatial Harmonics	Technique for parallel imaging based upon k-space calculations.
SORS-STC	Slice-selective Off-Resonance Sinc Pulse Saturation transfer contrast	Method for magnetization transfer imaging (Toshiba).
SPACE	Sampling Perfection with Application of optimized Contrasts using different flip angle Evolution	3D TSE sequence with variable refocusing flip angles (Siemens).
SPEEDER	SPEEDER	An image-based parallel imaging method related to SENSE (Toshiba; 'SPEEDER' is not an acronym, it is a tradename).
SPGR	SPoiled GRAdient echo	A spoiled GE imaging sequence using RF spoiling (GE Healthcare). FSPGR is a fast version of SPGR.
SPL	Sound Pressure Level	Acoustic noise, often measured with reference to the response of the human ear in decibels, dB(A).
SR	Slew Rate	The maximum gradient amplitude divided by its minimum rise time in $T\ m^{-1}\ s^{-1}$.

(cont.)

Acronym/symbol	Term	Explanation
SS	Slice Selection	A process using a slice-select gradient at the same time as a narrow-bandwidth RF pulse to excite a slice of tissue.
SS-TSE	Single Shot Turbo Spin Echo	A single shot TSE imaging sequence that also uses half-Fourier data acquisition and reconstruction. See also HASTE.
SSP	Sloped Slab Profile	Linear flip angle variation used in 3D-MRA to reduce saturation effects in slice-encoding direction (Hitachi).
STIR	Short TI Inversion Recovery	An IR-based imaging sequence used to null the signal from fat based upon its longitudinal relaxation time T_1 .
SWAN	T2-Star Weighted ANgiography	Method for susceptibility-weighted imaging (GE Healthcare).
SWI	Susceptibility-Weighted Imaging	Method for susceptibility-weighted imaging (Siemens).
SWIp	Susceptibility-Weighted Imaging with phase difference	Method for susceptibility-weighted imaging (Philips).
T_1	Spin-lattice relaxation time	Spin-lattice relaxation time, also known as longitudinal relaxation time. Characterizes the recovery of the longitudinal magnetization M_z towards M_0 .
T1-FFE		A spoiled GE imaging sequence using RF spoiling (Philips).
T_2	Spin-spin relaxation time	Spin-spin relaxation time, also known as transverse relaxation time. Characterizes the decay of transverse magnetization M_{xy} to zero.
T_2^*	Apparent Spin-spin relaxation time	Apparent spin-spin relaxation time. Characterizes the decay of the free induction signal to zero.
TA		Acquisition or scan time.
TE	Echo time	Time to the peak MR signal from the initial excitation.
TFE	Turbo Field Echo	A fast GE imaging sequence (Philips). See also FFE.
THRIVE	T_1 High Resolution Isotropic Volume Excitation	3D T_1 w GE with fat suppression used for dynamic liver imaging (Philips). The eTHRIVE variant (enhanced THRIVE) uses asymmetric k-space ordering. See also LAVA, VIBE
TI	Inversion Time	See IR (inversion recovery).
TIGRE	T_1 GRadient Echo with RF fat saturation	3D T_1 w GE with fat suppression used for dynamic breast and abdomen imaging (Hitachi).
Time-SLIP	Time-Spatial Labelling Inversion Pulse	Method for non-contrast MR angiography based on fully rewound GE (Toshiba).
TOF	Time of Flight	Non-contrast MRA technique that relies on the in-flow of unsaturated spins.

(cont.)

Acronym/symbol	Term	Explanation
TONE	Tilted Optimized Non-saturating Excitation	Linear flip angle variation used in 3D-MRA to reduce saturation effects in slice-encoding direction (Philips, Siemens). Also a generic term.
TR	Repetition time	Repetition time of an MR pulse sequence. Time between successive excitations.
TRANCE	TRiggered Angiography Non Contrast Enhanced	Method for non-contrast MR angiography based on 3D TSE (Philips).
TRAQ	Time-Resolved AcQuisition	Method for time-resolved CE-MRA (Hitachi).
TRICKS	Time-Resolved Imaging of Contrast KineticS	Method for time-resolved CE-MRA (GE Healthcare).
TRSG	Time-Reversed SARGE	A time-reversed GE imaging sequence (Hitachi).
TSE	Turbo Spin Echo	A multiple spin-echo sequence that uses a segmented k-space.
TWIST	Time-resolved angiography With Stochastic Trajectories	Method for time-resolved CE-MRA (Siemens).
UTE	Ultra-short TE	Ultra-short TE sequence capable of measuring the bound fraction of nuclei.
VASC ASL	VASCular Arterial Spin Labelling	Method for non-contrast MR angiography based on fully rewound GE (Hitachi).
VASC FSE	VASCular Fast Spin Echo	Method for non-contrast MR angiography based on 3D TSE (Hitachi).
VCG	Vector CardioGram	Method for ECG gating that helps reduce T-wave mistriggering.
v_{enc}		The maximum velocity that can be uniquely encoded in phase-contrast MR angiography.
VIBE	Volume Interpolated Breath-hold Examination	3D T ₁ w GE with fat suppression used for dynamic liver imaging (Siemens).
VIBRANT	Volume Image Breast Assessment	3D T ₁ w GE with fat suppression used for dynamic breast imaging (GE Healthcare).
VIEW	Volume Imaging with Echo Weighting	3D TSE sequence with variable refocusing flip angles (Philips).
WFS	Water–Fat Shift	The spatial misregistration between water and fat in the frequency encoding direction (Philips).

